




Captek provides one of the highest quality and most esthetic and biocompatible restorations available for your patients.

Advantages:

- Eliminates Gray Margins
- Extreme Bond Strength
- No special prep designs
- No Special Chair-side or bonding techniques required
- Superior Fit

Captek is an advanced metal-lurgic technology for producing high gold metal substructures for crowns and bridges. It is an internally reinforced gold material that exhibits both the color, workability and nobility of pure gold along with the strength and stability of the platinum/palladium family of cast alloys.



Today's patients are better informed and looking for quality dental health care. Captek's patented advanced gold material produces a crown that will look natural at the gum line, eliminating any unsightly gray or black lines.

**If you have any questions call
888-875-6011**