

Menu of Implant Abutment and Material Choices



		Straumann				Dentsply Friadent		BioHorizons®			Lifecore				Sybron Innova			BTI®		Thomen	
		Stock Titanium	Custom Cast	Procera®	Atlantis™	Stock Titanium	Custom Cast	Stock Titanium	Custom Cast	Atlantis™	Stock Titanium	Stock Zirconia	Atlantis™	UCLA	Stock Titanium	Custom Cast	Atlantis™	Stock Titanium	Custom Cast	Stock Titanium	Custom Cast
Connection	Internal	✓	✓	✓	✓	✓	✓													✓	✓
	Both							✓	✓		✓	✓		✓	✓			✓	✓		
	External								✓				✓			✓					
Abutment	Stock	✓				✓		✓			✓				✓			✓		✓	
	Custom Cast		✓				✓		✓				✓		✓				✓		✓
	CAD/CAM			✓	✓				✓			✓				✓					
Material	Gold Nitrate Titanium				✓					✓		✓				✓					
	Titanium	✓		✓	✓	✓		✓		✓		✓		✓		✓		✓		✓	
	Zirconia			✓						✓	✓	✓				✓					
	Alloy		✓				✓		✓				✓		✓				✓		✓
Final Screw	Titanium	✓	✓		✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
	Gold																	✓	✓		
	Torque Tight (Black)			✓																	