

Renew/refurbish your patient's comfortable old partial; add clasps and saddle areas. replace extracted teeth. Make old partials look like new and save chair time!

Case Study



8 & 9 were recently extracted. Patient liked the original partial. Teeth needed to be added to existing partial. Anterior saddle was directly cast and laser welded to partial, then acrylic and teeth added, saving the patient the expense of a new partial.

Expertec is a solder-free laboratory. Laser welding results in stronger and more accurate bridgework. It is cleaner, faster and stronger than conventional soldering techniques. Welding is more compatible in the oral environment.

- Partial denture reconstructions
- Implant and precision attachments
- Fabrication and repair of dental prosthetics
- Laser weld right on the master model ensures precise fit
- Repair / reassemble heat sensitive parts without the risk of warping
- Repair without disturbing acrylics, ceramics or porcelain materials
- Seam weld parts without the use of solder

Call for more information
888-875-6011