



DR. _____ DATE _____

TECHNICAL PREFERENCES

We customize your restorations based on your preferences. We have found most of our clients prefer their case fabricated in a similar fashion. To make communication a little easier we have developed a list of the most requested procedures. If you routinely prefer other options please check these and we will enter these preferences in our computer. When you send a prescription with your case any specific requests will take precedence over your routine preferences.

When sending in cases to us please:

1. Fill out laboratory scripts completely.
2. Send as much information as possible.
3. Include study models, diagnostic wax ups, temporary models, pictures and patient concerns.

COMMON PREFERENCES

- White High Noble 54%
- 3 Coats Die Spacer
- Narrow Lingual Band
- Porcelain to Buccal Margin
- Porcelain Occlusal
- Normal Embrasure
- Normal Proximal Contact
- Out of Occlusal .1MM
- Light Occlusal staining
- Normal Emergence Profile
- No Cervical Staining
- Saddle Tip Pontic (anterior)
- Partial Ridge (posteriors)
- Pontic Slight Cast Relief
- Slight Relief of Opposing Model for Clearance, if needed
- Contacts Double Checked on Solid Model

OPTIONS

- Yellow High Noble (86% gold) White Noble (gold & palladium)
- No spacer 1 coat 2 coats
- Other _____
- Narrow Band On Anterior
- On Posterior Porcelain Butt Margin
- Metal Occlusal On Anterior On Posterior
- Open Closed
- Broad Light Heavy
- .3MM .5MM Shimstock
- None Medium Heavy
- Convex Straight
- Light Heavy
- Bullet Full Ridge Ovoid Socket Partial Ridge
- No Relief Ovoid
- No Call Call first
- No Solid Model

May we routinely use a resin coping for slight prep reduction if needed? Yes No Call first

Please complete and send with your first case or fax to 734-641-4225